Ohaus School Balance
Model SB1200

- OPERATING MANUAL -

- Accurate and simple to use
- Durable lightweight plastic for easy clean up and long life

General Product Features

- **Interchangeable Red Pans:** Allows user to place a variety of objects for weighing. Removable for easy clean up.
- **Tare Poise:** Located on the rear of the beam, allows user to zero balance before weighing.
- **Recessed handles in base:** Allows the transport of the balance without disturbing the contents.

Balance Setup

- Remove balance from carton.
- Make sure that pans are placed securely on pan supports.
- Slide poise until pointer is lined up with center of graduated lines.

Weighing

To determine weight of object

- Make sure that pointer is lined up with center graduation.
- Place object to be weighed on one pan. The pointer will swing to the opposite side.
- Place masses on other pan until the pointer is again lined up with the center graduation.

To compare the weights of two objects

- Make sure that pointer is lined up with center graduation.
- Place first object to be weighed on one pan. The pointer will swing to the opposite side.
- Place second object on other pan. The pointer will point away from the heavier object.

Available Accessories

- 8 piece mass set 50g x 1g (80850117)
- 12 piece mass set 500g x 1g (80850114)
- 13 piece mass set 1000g x 1g (80850115)

EXPERIMENTS

Hot Air/Cool Air

Objective

To determine which is heavier – hot air or cool air.

Materials

- Two 8½' x 11" cardboard tablet backs
- Hair dryer or 75+ watt lamp bulb
- Small paper bags
- Ohaus SB1200 School Balance

Procedure

1. Place the cardboard pieces on each balance pan with a paper bag upside down on each piece of cardboard. Zero the balance.
2. Remove one of the bags and place over the lamp bulb for 2-3 minutes, being careful not to have the bulb touch the bag.

or

Using a hair dryer on high heat, fill the bag with hot air.
3. Without tipping the bag, place it back on the balance and observe the zero point.

Notes

- Which is heavier, hot air/cool air?
- Why is this so? (molecular motion)
- Talk about hot air balloons, ceiling fans, furnaces in the basement, molecular motion (solids, liquids, gasses)

Sugar in Soda

Objective

To show the relative amount of sugar in soda

Materials

- Same size can or bottle of regular and diet soda
- Sugar packets
- Ohaus SB1200 School Balance

Procedure

1. Place regular soda on the left pan and the diet soda on the right pan.
2. Add sugar packets to the right pan until the balance returns to a neutral position.
3. Remove the soda cans or bottles from the balance.
4. Empty sugar from the packets on to the left pan and weigh the sugar.

Note

The amount of error you will notice from your weighing and that which is on the can or bottle is caused by the following. The sugar in the soda is high fructose corn syrup. The diet contains aspartame and benzoic acid.

Limited Warranty

Ohaus products are warranted against defects in materials and workmanship for a period of 12 months from the date of delivery. During the warranty period Ohaus will repair, or at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.


The warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

Ohaus Corporation
7 Campus Drive
Suite 310
Parsippany, NJ 07054 USA
Tel: (973) 377-9000; Fax: (973) 944-7177

With offices worldwide
www.ohaus.com

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P/N 80250832 A © Ohaus Corporation 2010, alle Rechte vorbehalten.
P/N 80250832 A © Ohaus Corporation 2010, Tutti I diritti sono riservati.